

DynaVue⁺

Easy-to-use 2-in-1 Digital X-ray and Fluoroscopy System

Designed by JPI, DynaVue is an easy-to-use 2-in-1 digital radiography system designed for great X-ray images of companion animals. In development of DynaVue, our team worked side-by-side with veterinarians to optimize the system and streamline essential processes. DynaVue is also powered by ExamVue software, which was developed specifically with veterinarians and their workflow in mind.



Superior Digital Imaging and Increased Efficiency

- 17" x 17" Field of View
- 4-way floating tabletop for X-ray positioning with minimal disturbance to subject
- Easy-to-use operating console with large touchscreen installed on tube stand
- High sensitivity Cesium Iodide (CsI) detector technology
- Low dose, high frequency generator and remote capture feature
- Compact unit with movable casters for convenient positioning throughout your facility
- 220V, 15A power requirements for easy installation in any location.

Added Value of Fluoroscopy

You may not be utilizing fluoroscopy in your practice now, but there are plenty of reasons to expand your horizons. Improve your service and generate revenue performing procedures you couldn't do in the past. With the push of a button, this space-saving system can go from digital X-ray to fluoroscopy mode mid-exam, without any adjustment to the animal. Applications beyond digital X-ray are plentiful with DynaVue's added value fluoroscopy capabilities. Some of these procedures include:

- Minimally invasive fracture repair
- Intraluminal stenting for tracheal collapse treatment
- Treating urethral and colonic obstructive disease
- Vascular procedures
- Orthopedic procedures
- Contrast studies of the gastrointestinal tract



Contact Us:

(787) 462-6868 / (787) 645-1301
wmedsolutions@gmail.com